



## INJECTOR PORT LINERS - VARIAN GCS

All injection port liners are polysiloxane deactivated and ready for use. They are manufactured to exacting tolerances and come in various configurations for many popular GCs. Guaranteed inert for phenols, organic acids and bases.

CONFIGURATION		DIMENSIONS	OEM PART #	QTY	QUADREX CAT. #	PRICE
	Open Split 1075/1077	4mm I.D. X 72mm	16-000830-00	1	DGL-208311	\$24.26
				10	DGL-208311-10	\$195.45
	Frit Split 1075/1077	4mm I.D. X 72mm	01-900109-03 16-000830-01	1	DGL-208312	\$40.88
				10	DGL-208312-10	\$339.54
	Open Split w/ Glass Wool packed 1075/1077	4mm I.D. X 72mm	01-900109-01	1	DGL-208313	\$29.53
				10	DGL-208313-10	\$238.48
	Splitless Borosilicate 1075/1077	2mm I.D. X 74mm	01-900109-05 03-949437-00	1	DGL-208314	\$28.39
				10	DGL-208314-10	\$220.31
	SPI 1093/1094	0.5mm I.D. X 54mm	01-900109-06 03-918332-01	1	DGL-208316	\$42.01
				10	DGL-208316-10	\$310.01
	SPI 1093/1094	0.8mm I.D. X 54mm	01-900109-07 03-918332-02	1	DGL-208318	\$42.01
				10	DGL-208318-10	\$310.01